


Surgitex multiple choice Prostrate Prone Positioning

S1101-3

Foam cored gel
28 x 24 x 14 Cm



S1101-2

Basic gel (post operation)
25 x 20 x 10 Cm



The see-through hole

allows the pads to be used for wider applications. The canals provided on both sides allow an easy path for anaesthesia tubes


S1101-3

Two layer version
comprises a layer of
firmer gel base and
top layer of softer gel
25 x 20 x 11.5 Cm




S1101-4

Single layer soft gel
25 x 20 x 6 Cm



S1101-14

Single use foam with Acrylic detachable mirror 30 x 24 x 15 Cm



Surgitex Ltd.

96, Haven Chase, Cookridge, Leeds LS16 6SG
Tel: +44 8450 573288. Email: enquiries@surgitex.com

www.surgitex.com

